

# **FDA Submission**

**Your Name:** Sascha Metzger

**Name of your Device:** Chest X-Rays Pneumonia Detector

## **Algorithm Description**

### **1. General Information**

**Intended Use Statement:** Help Radiologists in detecting Pneumonia in Chest X-Rays

**Indications for Use:**

- Applicable for men and women from 1 to 90 years old
- Chest X-Ray image must be taken in the AP or PA position
- Chest X-Ray image must be in DICOM format


**Device Limitations:**

- Diagnosis can be made on a computer with a standard CPU, although a GPU is preferred

**Clinical Impact of Performance:**

- The model has a lower precision and higher recall
- This means the model is most confident when the test result is negative
- Therefore it is best used for worklist prioritization and not as a diagnosis tool

### **2. Algorithm Design and Function**



#### **DICOM Checking Steps:**

- Modality is "DX"
- Body part examined is "CHEST"
- Patient Position is "PA" or "AP"

#### **Preprocessing Steps:**

- Image is normalized
- Image is reshaped
- Image is repeated across 3 channels

#### **CNN Architecture:**

Model: "sequential\_4"

Layer (type)	Output Shape	Param #
model_4 (Model)	(None, 7, 7, 512)	14714688
flatten_4 (Flatten)	(None, 25088)	0
dropout_11 (Dropout)	(None, 25088)	0
dense_11 (Dense)	(None, 1024)	25691136
dropout_12 (Dropout)	(None, 1024)	0
dense_12 (Dense)	(None, 512)	524800
dropout_13 (Dropout)	(None, 512)	0
dense_13 (Dense)	(None, 256)	131328
dropout_14 (Dropout)	(None, 256)	0
dense_14 (Dense)	(None, 1)	257
Total params: 41,062,209		
Trainable params: 28,707,329		
Non-trainable params: 12,354,880		


- The model is based on the VGG16 model
- The model uses the first 17 layers of the VGG16 model
- The VGG16 model output is flattened and passed through several additional dense and dropout layers

### 3. Algorithm Training

#### Parameters:

- Types of augmentation used during training
  - Horizontal flips
  - No Vertical flip s
  - Height shift range of 0.1,
  - Width shift range of 0.1,
  - Rotation range of 25,
  - Shear range of 0.1,
  - Zoom range of 0.15
- Batch size: 64
- Optimizer learning rate: 1e-4
- Layers of pre-existing architecture that were frozen: First 17 layers of VGG model

- Layers of pre-existing architecture that were fine-tuned: None
- Layers added to pre-existing architecture: Flatten, Dense and Dropout layers



#### Final Threshold and Explanation:

- Threshold: 0.5111328
- F1 Score: 0.3713298791018998

The final threshold of 0.51 was based on the highest F1 Score of 0.37. Based on [this paper](#) the average radiologist has a F1 Score of 0.387. This means this model achieved nearly the same performance as the average radiologist.

#### **4. Databases**

- The Dataset can be found here: [NIH Chest X-ray Dataset](#)
- It contains 112,120 chest x-ray images
- Each image has the following meta data:
  - Image Index
  - Finding Labels
  - Follow-up #
  - Patient ID
  - Patient Age
  - Patient Gender
  - View Position
  - Original Image Size
  - Original Image Pixel Spacing

#### **Description of Training Dataset:**

- The training data is split equally between Pneumonia and non Pneumonia patients
- It contains 2290 images

#### **Description of Validation Dataset:**

- The training data has 20% Pneumonia and 80% non Pneumonia patients
- It contains 1430 images

#### **5. Ground Truth**

- 112,120 X-ray images with disease labels from 30,805 unique patients
- The disease labels were created using Natural Language Processing (NLP) to mine the associated radiological reports
- The labels are expected to be >90% accurate and suitable for weakly-supervised learning
- The data includes 14 common thoracic pathologies:
  - Atelectasis
  - Consolidation
  - Infiltration
  - Pneumothorax

- Edema
- Emphysema
- Fibrosis
- Effusion
- Pneumonia
- Pleural thickening
- Cardiomegaly
- Nodule
- Mass
- Hernia

## **6. FDA Validation Plan**

### **Patient Population Description for FDA Validation Dataset:**

- Applicable for men and women from 1 to 90 years old
- Chest X-Ray image must be taken in the AP or PA position
- Chest X-Ray image must be in DICOM format

### **Ground Truth Acquisition Methodology:**

- Silver Standard: Validation by 3 different radiologists

### **Algorithm Performance Standard:**

- The algorithm's F1 score should be more than that of average radiologist (0.387)